

Compare the i2 to the i3 - Quick Glance



Replaced when battery runs out (3-5 years)

Sensitive to cell phone signals

"Can we trust the data?"

Always on

Connect to computer using the dosimeter cradle

Non-replaceable metal pocket clip

Colored paper inlay rings, no obvious space for name

Battery replaceable – increased life cycle (<10 years)

Better EMC immunity

Improved measurement specification

Automatic sleep mode

Direct USB cable connection to computer

User replaceable pocket clip, made from reinforced plastic

Colored paper inlays with space for name

