



#### PRODUCT FEATURES:

- Customized for use with complementary RaySafe Pro-Slit Camera and the X-ray Pinhole Assemblies

## RaySafe Pro-Stand

Adjustable stand, designed to make focal spot measuring procedures easy to perform, as well as ensuring accurate results. It can be used with our pinholes and slit camera (not included) and can also be used for manual HVL measurements.

#### TECHNICAL DATA

- Adjustable height from 350-600 mm (11.8-23.6 in.), providing a wide magnification range
- Adjustable horizontal position and vertical angle of the table top
- Optional positioning tool for easy setup verification – especially useful for measurements on mammography units

#### BASIC VERSION

Product No. 4829231

- Adjustable height
- Adjustable horizontal position
- Compatible with our pinholes and slit camera
- Comes in cardboard with fitted inlay

#### FULL VERSION

Product No. 4829317

- All features from Basic version
- Possibility to tilt with an angle of 10 degrees for easier positioning on mammography systems
- Includes positioning tool for easy setup
- Comes in a heavy duty rugged case for safe transportation